

# Antibody ELISA Assays

## Monitor 28 Common Pathogens

- Mouse
- Rat
- Rabbit
- Guinea Pig
- Hamster

## Ease of Use

- Standard ELISA format
- Easy to follow directions, no special training needed
- Common protocol for all assays
- Compatible with plate washers and plate readers
- Ready-to-use ELISA plates, conjugate and substrate

## Reliability


- Reproducible: CV <5%
- Detects low levels of antibody for early identification of infection
- Early, reliable test results to insure good laboratory animal health

## Cost-Effectiveness

- Significant savings vs. outside testing service
- 2 hour turnaround time
- 12 month shelf life
- Unused strips can be removed and used later
- Components available separately
- Shipping and Storage at 4°C

## Multiplexing

- Up to six different pathogens per well
- Custom made plates to fit your need
- Drastically reduce labor and time



## PCR Assays

- PCR Controls for 28 Common Pathogens
- Helicobacter hepaticus PCR Kit

## DNA Extraction Kits

- Mouse Tail
- Mouse Hair
- Eartags



503 Gateway Drive West  
Thurmont, MD 21788 USA  
Phone: +301-228-2444  
Fax: +301-560-6570  
US Toll Free: 888-562-8914

[www.xpressbio.com](http://www.xpressbio.com)

All Products and Services for Research Use Only