

LAB ANIMAL HEALTH DIAGNOSTICS

Monitoring antibodies to Common Lab Animal Pathogens

Multiplex & Single Agent ELISAs

- Up to six different pathogens per well
- Customized to meet specific testing needs
- Significant reduction in labor and turnaround time
- 5ul sample per assay

Easy to Use

- 96 well ELISA plates
- One protocol for all assays
- Ambient Shipping, Storage at 2-8°C
- Compatible with plate washers, readers & automated plate processors
- Ready-to-use ELISA plates, controls, diluent, conjugate and substrate



Reliable

- Highly Reproducible Results: CVs < 5%
- Low levels of antibody detected for early identification of infection
- Consistent performance in Stability Testing at 4°C and 37°C

Cost-Effective

- 18 month shelf life
- 2 and ½ hour turnaround time
- Unused strips can be removed and used later
- Components available separately and in bulk
- Shipping at ambient temperature (no dry ice)

PCR Assays

- PCR Controls and Assays
- Helicobacter genus & 6 species

DNA Extraction Kits

- Mouse Tail
- Ear Tags
- Mouse Hair



503 Gateway Drive West
Thurmont, MD 21788 USA
Phone: 301-228-2444
Fax: 301-560-6570
US Toll Free: 888-562-8914

WWW.XPRESSBIO.COM