

HIV PRODUCTS

HIV-1 Integrase Assay

- Measures effects of interacting proteins
- Measures effects of inhibition of HIV-1 integrase activity
- Quantitative
- Simple same day assay

HIV-1 Integrase Protein

- Full length recombinant
- Purified from bacteria

HIV-1 p24 ELISA Assay

- Tissue culture supernatant
- Linear up to 200pg/ml
- Standard ELISA format
- Easy to follow directions, no special training needed
- Easily automatable

HIV and SIV Antibodies and Recombinant Proteins

- Monoclonals
- Polyclonals
- Recombinants
- *E. coli* or Baculovirus expressed
- Biotin, FITC, HRP labels
- Frozen or lyophilized



503 Gateway Drive West
Thurmont, MD 21788 USA
Phone: +301-228-2444
Fax: +301-560-6570
US Toll Free: 888-562-8914

www.xpressbio.com